


Anti-Thop1 Rabbit Polyclonal Antibody For Wb, Elisa - P52888

Item Number: CM0030867



Introduction

Rabbit polyclonal antibody against human THOP1, validated for WB and ELISA; reacts with human, mouse, rat; detects 79kDa thimet oligopeptidase involved in neuropeptide metabolism and amyloid-beta degradation.

[Learn More](#)

Parameter	Value
Product Name	THOP1 Rabbit pAb
Remarks/Alias	TOP; MP78; EP24.15; MEPD_HUMAN; THOP1
Species	Human
Gene ID (Human)	7064
Gene ID	7064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human THOP1 (NP_003240.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52888
Protein Weight	79kDa
Shipping	Ice bag