

Anti-Thra Polyclonal Antibody For Wb, Ihc-P, Elisa - P10827

Item Number: CM0030230



Introduction

Rabbit polyclonal antibody against THRA (Thyroid hormone receptor alpha), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for nuclear receptor transcription studies.

[Learn More](#)

Parameter	Value
Product Name	THRA Rabbit pAb
Remarks/Alias	AR7; EAR7; ERBA; CHNG6; ERBA1; NR1A1; THRA1; THRA2; c-erbA; ERB-T-1; TRalpha; THRalpha; TRalpha1; TRalpha2; c-ERBA-1; THRalpha1; THRalpha2; THRA
Species	Human
GeneID (Human)	7067
GeneID	7067
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10827
Protein Weight	55kDa
Shipping	Ice bag