

Anti-Thrap3 Polyclonal Antibody For Wb - Q9Y2W1

Item Number: CM0029955

Introduction

Rabbit polyclonal antibody against THRAP3 (TRAP150), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying pre-mRNA splicing, mRNA decay, and circadian clock regulation.

[Learn More](#)

Parameter	Value
Product Name	THRAP3 Rabbit pAb
Remarks/Alias	BCLAF2, TRAP150, THRAP3
Species	Human
Gene ID (Human)	9967
Gene ID	9967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-500 of human THRAP3 (NP_005110.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y2W1
Protein Weight	109kDa
Shipping	Ice bag