

Anti-Thromboxane Synthase Rabbit Monoclonal Antibody For Wb And Elisa - P24557

Item Number: CM0008345

Introduction

Rabbit monoclonal antibody to human Thromboxane synthase (TXAS/CYP5A1) validated for WB and ELISA. Cross-reacts with mouse. Ideal for thromboxane signaling and platelet function studies.

[Learn More](#)

Parameter	Value
Product Name	Thromboxane synthase Rabbit mAb
Remarks/Alias	TS; TXS; CYP5; THAS; TXAS; CYP5A1; GHOSAL; BDPLT14; Thromboxane synthase
Species	Human
GeneID (Human)	6916
GeneID	6916
Immunogen	Synthetic Peptide / A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Thromboxane synthase (P24557)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24557
Protein Weight	61kDa
Shipping	Ice bag