

Anti-Thromboxane Synthase Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24557

Item Number: CM0023216

Introduction

Thromboxane synthase Rabbit Polyclonal Antibody (P24557) validated for WB, IHC-P, ELISA in human, mouse, rat. Suitable for studying thromboxane A2 synthesis, platelet aggregation, and vasoconstriction.

[Learn More](#)

Parameter	Value
Product Name	Thromboxane synthase Rabbit pAb
Remarks/Alias	TS; TXS; CYP5; THAS; TXAS; CYP5A1; GHOSAL; BDPLT14; Thromboxane synthase
Species	Human
GeneID (Human)	6916
GeneID	6916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 355-534 of human Thromboxane synthase (NP_001052.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24557
Protein Weight	61kDa
Shipping	Ice bag