

Anti-Thrsp Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q92748

Item Number: CM0029423

Introduction

THRSP Rabbit Polyclonal Antibody for WB, IF/ICC, ELISA detects thyroid hormone-inducible hepatic protein. Cross-reacts with human, mouse, rat. Essential tool for lipogenesis and transcription research. Q92748.

[Learn More](#)

Parameter	Value
Product Name	THRSP Rabbit pAb
Remarks/Alias	S14; Lpgp; THRP; LPGP1; SPOT14; THRSP
Species	Human
Gene ID (Human)	7069
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-146 of human THRSP (NP_003242.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92748
Protein Weight	17kDa
Shipping	Ice bag