


Anti-Thymidine Kinase 1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P04183

Item Number: CM0026729



Introduction

Anti-Thymidine Kinase 1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA with human, mouse, rat reactivity. Thymidine kinase 1 is a cell cycle-regulated enzyme essential for nucleotide metabolism and prodrug activation. Suitable for cancer and cell proliferation research.

[Learn More](#)

Parameter	Value
Product Name	Thymidine Kinase 1 Rabbit pAb
Remarks/Alias	TK2; Thymidine Kinase 1
Species	Human
Gene ID (Human)	7083
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-194 of human Thymidine Kinase 1 (NP_003249.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04183
Protein Weight	25kDa
Shipping	Ice bag