

Anti-Thymidylate Synthase (Tsy) Polyclonal Antibody, Ko Validated - P04818

Item Number: CM0024110



Introduction

CamelBio offers high-quality rabbit polyclonal antibody against human thymidylate synthase (TYSY). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for nucleotide metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TYSY Rabbit pAb
Remarks/Alias	TS; TMS; DKCD; HST422; SY
Species	Human
GeneID (Human)	7298
GeneID	7298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-313 of human TYSY (NP_001062.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04818
Protein Weight	36kDa
Shipping	Ice bag