


Anti-Thyroglobulin Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P01266

Item Number: CM0005071



Introduction

This rabbit monoclonal antibody recognizes human Thyroglobulin (~305 kDa) and is validated for immunofluorescence (IF-P), immunohistochemistry (IHC-P), and ELISA. It cross-reacts with mouse and rat, making it suitable for thyroid hormone biosynthesis studies and diagnostic test development.

[Learn More](#)

Parameter	Value
Product Name	Thyroglobulin Rabbit mAb
Remarks/Alias	TGN; AITD3; Thyroglobulin
Species	Human
GeneID (Human)	7038
GeneID	7038
Immunogen	Synthetic peptide corresponding to amino acids 1000-1100 of human Thyroglobulin (P01266)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01266
Protein Weight	305kDa
Shipping	Ice bag