

Anti-Thyroid Stimulating Hormone (Tsh) Monoclonal Antibody For If, Ihc-P, Elisa - P01215/P01222

Item Number: CM0005430



Introduction

Rabbit monoclonal antibody targeting human Thyroid Stimulating Hormone (TSH) alpha and beta subunits. Suitable for IF, IHC-P, and ELISA applications. Ideal for research on hormone signaling and pituitary function.

[Learn More](#)

Parameter	Value
Product Name	Thyroid Stimulating Hormone Rabbit mAb
Remarks/Alias	HCG; LHA; FSHA; GPA1; GPHa; TSHA; GPHA1; CG-ALPHA; TSHB
Species	Human
GeneID (Human)	1081/7252
GeneID	1081/7252
Immunogen	Recombinant protein corresponding to amino acids 27-116 of human CGA and 21-138 of human TSHB (NP_000726.1/NP_000540.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01215/P01222
Protein Weight	13kDa/16kDa
Shipping	Ice bag