

Anti-Thyroid Stimulating Hormone (Tsh) Polyclonal Antibody For Ihc-P, Elisa - P01215 / P01222

Item Number: CM0020564



Introduction

Rabbit polyclonal antibody against human Thyroid Stimulating Hormone (TSH), recognizing both alpha and beta subunits. Validated for IHC-P and ELISA, with cross-reactivity to mouse. Ideal for endocrine research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Thyroid Stimulating Hormone (TSH) Rabbit pAb
Remarks/Alias	HCG; LHA; FSHA; GPA1; GPHa; TSHA; GPHA1; CG-ALPHA; TSH-B; TSH-BETA; Thyroid Stimulating Hormone (TSH)
Species	Human
GeneID (Human)	1081/7252
GeneID	1081/7252
Immunogen	Recombinant protein corresponding to a sequence within amino acids 25-116&21-138 of human Thyroid Stimulating Hormone (TSH)(NP_000726.1/NP_000540.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	P01215/P01222
Protein Weight	13KDa/15KDa
Shipping	Ice bag