


Anti-Tia1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P31483

Item Number: CM0014057



Introduction

High-quality rabbit polyclonal antibody targeting human TIA1 (P31483). Suitable for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects human and mouse TIA1. Useful for alternative splicing, stress granules, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	TIA1 Rabbit pAb
Remarks/Alias	WDM; ALS26; TIA-1; TIA1
Species	Human
GeneID (human)	7072
GeneID	7072
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 287-386 of human TIA1 (NP_071505.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31483
Protein Weight	43kDa
Shipping	Ice bag