


Anti-Tial1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q01085

Item Number: CM0025971



Introduction

TIAL1 Rabbit pAb detects nucleolysin TIAR in human, mouse, and rat. Validated for western blot, immunohistochemistry, and ELISA. Suitable for alternative splicing and stress granule research.

[Learn More](#)

Parameter	Value
Product Name	TIAL1 Rabbit pAb
Remarks/Alias	TCBP; TIAR; TIAL1
Species	Human
Gene ID (Human)	7073
Gene ID	7073
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 266-375 of human TIAL1 (NP_003243.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01085
Protein Weight	42kDa
Shipping	Ice bag