


Anti-Tifa Monoclonal Antibody For Wb, Ihc-P, Elisa - Q793I8

Item Number: CM0011976



Introduction

KO validated anti-TIFA rabbit monoclonal antibody for WB, IHC-P, ELISA. Detects mouse TIFA, a 21 kDa adapter protein in NF-κB innate immune signaling. Ideal for infection and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TIFA Rabbit mAb
Remarks/Alias	T2bp
Species	Mouse
GeneID (Human)	92610
GeneID (Mouse)	211550
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-184 of mouse TIFA (NP_001343359.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse
SWISS	Q793I8
Protein Weight	21 kDa
Shipping	Ice bag