

Anti-Tifab Rabbit Monoclonal Antibody (Kd Validated) For Wb And Elisa - Q6Znk6

Item Number: CM0005554

Introduction

KD validated rabbit monoclonal antibody targeting human TIFAB, suitable for WB and ELISA applications. Ideal for investigating TRAF6 signaling and NF-kappaB pathway regulation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TIFAB Rabbit mAb
Species	Human
Gene ID (Human)	497189
Gene ID	497189
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 62-161 of human TIFAB(NP_001092691.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6ZNK6
Protein Weight	18KDa
Shipping	Ice bag