

Anti-Tify11B Polyclonal Antibody For Wb And Elisa - Q9C9E3

Item Number: CM0023558

Introduction

Arabidopsis thaliana TIFY11B (JAZ6) rabbit polyclonal antibody, targeting a jasmonate signaling repressor. Validated for Western blotting and ELISA. Ideal for plant defense and jasmonic acid pathway research.

[Learn More](#)

Parameter	Value
Product Name	TIFY11B Rabbit pAb
Remarks/Alias	jasmonate-zim-domain protein 6; T10D10.8; T10D10_8; TIFY DOMAIN PROTEIN 11B; TIFY11B
Species	Arabidopsis thaliana
GeneID	843577
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 170-269 of Arabidopsis thaliana TIFY11B (NP_565043.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q9C9E3
Protein Weight	30kDa
Shipping	Ice bag