

Anti-Tigar Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nq88

Item Number: CM0020157



Introduction

Rabbit polyclonal antibody against TIGAR (TP53-induced glycolysis and apoptosis regulator). Validated for WB, IF/ICC, ELISA in human, mouse, rat samples. Regulates glycolysis and apoptosis; used in cancer, stress, and metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TIGAR Rabbit pAb
Remarks/Alias	FR2BP; C12orf5; TIGAR
Species	Human
Gene ID (Human)	57103
Gene ID	57103
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human TIGAR (NP_065108.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ88
Protein Weight	30kDa
Shipping	Ice bag