


Anti-Tigit Rabbit Polyclonal Antibody For Wb, Elisa - Q495A1

Item Number: CM0022571



Introduction

High-quality rabbit polyclonal antibody targeting human TIGIT (VSIG9, VSTM3), validated for Western blot and ELISA applications. Cross-reacts with mouse TIGIT. Ideal for immune checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	TIGIT Rabbit pAb
Remarks/Alias	VSIG9; VSTM3; WUCAM; TIGIT
Species	Human
GeneID	201633 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 167-244 of human TIGIT (NP_776160.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q495A1
Protein Weight	26kDa
Shipping	Ice bag