


Anti-Tim-1/Kim-1/Havcr Rabbit Polyclonal Antibody - Q5Qns5

Item Number: CM0023656



Introduction

Rabbit polyclonal antibody targeting mouse TIM-1/KIM-1/HAVCR, suitable for WB and ELISA. Shows human cross-reactivity. Detects a 33 kDa protein, ideal for research in kidney injury, T-cell trafficking, and efferocytosis.

[Learn More](#)

Parameter	Value
Product Name	TIM-1/KIM-1/HAVCR Rabbit pAb
Remarks/Alias	Tim1; KIM-1; TIM-1; Timd1; TIM-1/KIM-1/HAVCR
Species	Mouse
Gene ID	171283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-237 of mouse TIM-1/KIM-1/HAVCR (NP_599009.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5QNS5
Protein Weight	33kDa
Shipping	Ice bag