


Anti-Tim-3/Havcr2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q8Tdq0

Item Number: CM0009934



Introduction

Rabbit monoclonal antibody against human TIM-3/HAVCR2. Validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse. Suitable for immune checkpoint and T-cell research. 33 kDa.

[Learn More](#)

Parameter	Value
Product Name	Anti-TIM-3/HAVCR2 Rabbit Monoclonal Antibody
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2
Species	Human
GeneID (Human)	84868
GeneID	84868
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag