

Anti-Tim-3/Havcr2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Tdq0

Item Number: CM0020107



Introduction

Rabbit polyclonal antibody against human TIM-3/HAVCR2, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for immune checkpoint research and T-cell exhaustion studies.

[Learn More](#)

Parameter	Value
Product Name	TIM-3/HAVCR2 Rabbit pAb
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2; TIM-3/HAVCR2
Species	Human
GeneID (Human)	84868
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-202 of human TIM-3/HAVCR2 (NP_116171.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag