

Anti-Tim23 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - O14925

Item Number: CM0008874

Introduction

Validated rabbit monoclonal antibody against TIM23 (TIMM23) for WB, IP, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Research use in mitochondrial import and mitophagy. 22kDa.

[Learn More](#)

Parameter	Value
Product Name	TIM23 Rabbit mAb
Remarks/Alias	TIM23
Species	Human
GeneID (Human)	100287932
GeneID	100287932
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 7-72 of human TIMM23 (NP_006318.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14925
Protein Weight	22kDa
Shipping	Ice bag