

Anti-Tim44 Rabbit Monoclonal Antibody For Wb And Elisa - O43615

Item Number: CM0003420



Introduction

Rabbit monoclonal antibody against human TIM44, a mitochondrial import inner membrane translocase subunit. Validated for WB, ELISA; cross-reacts with human, rat. Suitable for mitochondrial protein import research.

[Learn More](#)

Parameter	Value
Product Name	TIM44 Rabbit mAb
Remarks/Alias	TIM44
Species	Human
GeneID (Human)	10469
GeneID	10469
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-240 of human TIM44 (NP_006342.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O43615
Protein Weight	51kDa
Shipping	Ice bag