


Anti-Timeless Monoclonal Antibody For Wb, Elisa - Q9Uns1

Item Number: CM0002663



Introduction

Rabbit monoclonal antibody against human Timeless (TIM, hTIM) for Western blot and ELISA applications. Suitable for studying DNA replication, repair, and circadian clock regulation. Specific to human, molecular weight 139 kDa.

[Learn More](#)

Parameter	Value
Product Name	Timeless Rabbit mAb
Remarks/Alias	TIM; TIM1; hTIM; FASPS4; Timeless
Species	Human
GeneID (Human)	8914
GeneID	8914
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1058-1208 of human Timeless (Q9UNS1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UNS1
Protein Weight	139kDa
Shipping	Ice bag