


Anti-Timm10 Polyclonal Antibody For Wb, If/Icc, Elisa - P62072

Item Number: CM0029581



Introduction

Rabbit polyclonal antibody against human TIMM10 (mitochondrial import inner membrane translocase subunit Tim10), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial protein import studies.

[Learn More](#)

Parameter	Value
Product Name	TIMM10 Rabbit pAb
Remarks/Alias	TIM10; TIM10A; TIMM10A; TIMM10
Species	Human
Gene ID (Human)	26519
Gene ID	26519
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62072
Protein Weight	10kDa
Shipping	Ice bag