

Anti-Timm17B Monoclonal Antibody For Wb - O60830

Item Number: CM0007723

Introduction

TIMM17B Rabbit mAb is a monoclonal antibody targeting human TIMM17B, validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for mitochondrial import and translocation studies.

[Learn More](#)

Parameter	Value
Product Name	TIMM17B Rabbit mAb
Remarks/Alias	JM3; TIM17B; DXS9822; TIMM17B
Species	Human
GeneID (Human)	10245
GeneID	10245
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 60-160 of human TIMM17B (NP_005825.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60830
Protein Weight	18kDa
Shipping	Ice bag