

## Anti-Timm17B Polyclonal Antibody For Wb, Elisa - O60830

Item Number: CM0024316

### Introduction

High-quality TIMM17B rabbit polyclonal antibody validated for WB and ELISA in human and rat. Targets mitochondrial import inner membrane translocase subunit Tim17-B, essential for protein translocation across the mitochondrial inner membrane.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TIMM17B Rabbit pAb                                                                                                           |
| Remarks/Alias    | JM3; TIM17B; DXS9822; TIMM17B                                                                                                |
| Species          | Human                                                                                                                        |
| GeneID (Human)   | 10245                                                                                                                        |
| GeneID           | 10245                                                                                                                        |
| Immunogen        | Synthetic peptide - A synthetic peptide corresponding to a sequence within amino acids 53-152 of human TIMM17B(NP_005825.1). |
| Source           | Rabbit                                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                                        |
| Application      | WB, ELISA                                                                                                                    |
| Cross Reactivity | Human, Rat                                                                                                                   |
| SWISS            | O60830                                                                                                                       |
| Protein Weight   | 18kDa                                                                                                                        |
| Shipping         | Ice bag                                                                                                                      |