

Anti-Timm17B Rabbit Polyclonal Antibody For Wb, Elisa - O60830

Item Number: CM0028683

Introduction

TIMM17B Rabbit pAb validated for WB, ELISA. Detects human, mouse TIMM17B (18kDa), a mitochondrial inner membrane translocase. Ideal for mitochondrial import studies.

[Learn More](#)

Parameter	Value
Product Name	TIMM17B Rabbit pAb
Remarks/Alias	JM3, TIM17B, DXS9822, TIMM17B
Species	Human
GeneID (Human)	10245
GeneID	10245
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human TIMM17B (NP_005825.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O60830
Protein Weight	18kDa
Shipping	Ice bag