

Anti-Timm22 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Cq85

Item Number: CM0022298

Introduction

Anti-Timm22 rabbit polyclonal antibody validated for Western blot, immunoprecipitation, and ELISA. Reacts with human, mouse, rat Timm22, a mitochondrial inner membrane translocase subunit essential for multipass transmembrane protein import. Ideal for IVD research and development.

[Learn More](#)

Parameter	Value
Product Name	Timm22 Rabbit pAb
Remarks/Alias	Tim22; Timm22
Species	Mouse
Gene ID	56322
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 19-121 of mouse Timm22 (NP_062792.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9CQ85
Protein Weight	20kDa
Shipping	Ice bag