

## Anti-Timm23 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O14925

Item Number: CM0029387

### Introduction

Rabbit polyclonal antibody targeting human TIMM23 (TIM23), essential for mitochondrial protein import and mitophagy regulation. Validated for WB, IHC-P, IF/ICC, ELISA with cross-reactivity to mouse and rat.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | TIMM23 Rabbit pAb                                                                                                |
| Remarks/Alias    | TIM23; TIMM23                                                                                                    |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 100287932                                                                                                        |
| Gene ID          | 100287932                                                                                                        |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 7-72 of human TIMM23 (NP_006318.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | O14925                                                                                                           |
| Protein Weight   | 22kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |