

Anti-Timm29 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q8Bgx2

Item Number: CM0020204

Introduction

Rabbit polyclonal antibody targeting mouse Timm29 (SWISS Q8BGX2), suitable for WB, IP, ELISA. Cross-reacts with human. For mitochondrial protein import and TIM22 complex studies.

[Learn More](#)

Parameter	Value
Product Name	Timm29 Rabbit pAb
Remarks/Alias	TIM29; 1810026J23Rik; Timm29
Species	Mouse
GeneID	69773
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 157-266 of mouse Timm29 (NP_848734.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8BGX2
Protein Weight	29kDa
Shipping	Ice bag