

Anti-Timm44 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43615

Item Number: CM0024024



Introduction

TIMM44 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, ELISA with reactivity to human, mouse, rat. It targets TIMM44 (TIM44), a mitochondrial inner membrane translocase essential for protein import.

[Learn More](#)

Parameter	Value
Product Name	TIMM44 Rabbit pAb
Remarks/Alias	TIM44; TIMM44
Species	Human
Gene ID (Human)	10469
Gene ID	10469
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 153-452 of human TIMM44 (NP_006342.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43615
Protein Weight	51kDa
Shipping	Ice bag