

Anti-Timm8A Polyclonal Antibody For Wb, If/Icc, Elisa - O60220

Item Number: CM0028699

Introduction

Rabbit polyclonal antibody targeting human mitochondrial import inner membrane translocase subunit Tim8 A. Validated for Western Blot, IF/ICC, and ELISA; cross-reacts with mouse and rat. For mitochondrial transport and chaperone research.

[Learn More](#)

Parameter	Details
Product Name	TIMM8A Rabbit pAb
Remarks/Alias	DDP; MTS; DDP1; DFN1; TIM8; TIMM8A
Species	Human
GeneID (Human)	1678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of human TIMM8A (NP_004076.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60220
Protein Weight	11kDa
Shipping	Ice bag