

Anti-Timm8B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y5J9

Item Number: CM0024447

Introduction

Rabbit polyclonal antibody to TIMM8B (mitochondrial import inner membrane translocase subunit Tim8 B) for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for mitochondrial protein import and chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	TIMM8B Rabbit pAb
Remarks/Alias	DDP2; TIM8B; TIMM8B
Species	Human
Gene ID (Human)	26521
Gene ID	26521
Immunogen	Recombinant Protein, Recombinant fusion protein containing a sequence corresponding to amino acids 16-98 of human TIMM8B (NP_036591.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5J9
Protein Weight	9kDa
Shipping	Ice bag