

Anti-Timp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P01033

Item Number: CM0017433



Introduction

High-quality TIMP1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat TIMP1. Suitable for studying MMP inhibition, erythropoiesis, and integrin signaling in cancer and inflammatory diseases.

[Learn More](#)

Parameter	Value
Product Name	TIMP1 Rabbit pAb
Remarks/Alias	EPA; EPO; HCI; CLGI; TIMP; TIMP-1; TIMP1
Species	Human
GeneID (Human)	7076
GeneID	7076
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 128-207 of human TIMP1 (NP_003245.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01033
Protein Weight	23kDa
Shipping	Ice bag