

Anti-Timp2 Monoclonal Antibody For Wb, Elisa - P16035

Item Number: CM0003630



Introduction

Rabbit monoclonal antibody against human TIMP2 (Metalloproteinase inhibitor 2), validated for WB and ELISA. Cross-reacts with mouse, rat, and African green monkey. Suitable for MMP research.

[Learn More](#)

Parameter	Value
Product Name	TIMP2 Rabbit mAb
Remarks/Alias	DDC8; CSC-21K; TIMP2
Species	Human
Gene ID (Human)	7077
Gene ID	7077
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-220 of human TIMP2. (NP_003246.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, African green monkey
SWISS	P16035
Protein Weight	24kDa
Shipping	Ice bag