

Anti-Timp2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16035

Item Number: CM0013001



Introduction

Rabbit polyclonal antibody against human TIMP2, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying MMP inhibition, tissue remodeling, and cancer metastasis.

[Learn More](#)

Parameter	Value
Product Name	TIMP2 Rabbit pAb
Remarks/Alias	DDC8; CSC-21K; TIMP2
Species	Human
GeneID (human)	7077
GeneID	7077
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human TIMP2 (NP_003237.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16035
Protein Weight	24kDa
Shipping	Ice bag