

Anti-Timp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P16035

Item Number: CM0021530



Introduction

Anti-TIMP2 rabbit polyclonal antibody for detecting metalloproteinase inhibitor 2. Tested in WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for ECM and metastasis research.

[Learn More](#)

Parameter	Value
Product Name	TIMP2 Rabbit pAb
Remarks/Alias	DDC8; CSC-21K; TIMP2
Species	Human
GeneID (human)	7077
GeneID	7077
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human TIMP2 (NP_003246.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16035
Protein Weight	24kDa
Shipping	Ice bag