

Anti-Timp3 Rabbit Polyclonal Antibody For Wb And Elisa - P35625

Item Number: CM0013092



Introduction

Rabbit polyclonal antibody against TIMP3, validated for WB and ELISA. Detects human, mouse, and rat TIMP3 (24 kDa). Ideal for extracellular matrix, cancer, and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	TIMP3 Rabbit pAb
Remarks/Alias	SFD; K222; K222TA2; HSMRK222; TIMP3
Species	Human
GeneID (Human)	7078
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 110-210 of human TIMP3 (NP_000353.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35625
Protein Weight	24kDa
Shipping	Ice bag