

Anti-Timp4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99727

Item Number: CM0026819

Introduction

TIMP4 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for extracellular matrix remodeling, cardiovascular biology, and MMP inhibition studies.

[Learn More](#)

Parameter	Value
Product Name	TIMP4 Rabbit pAb
Remarks/Alias	TIMP-4; TIMP4
Species	Human
GeneID	7079
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-224 of human TIMP4 (NP_003247.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99727
Protein Weight	26kDa
Shipping	Ice bag