

Anti-Tin2/Tinf2 Polyclonal Antibody For Wb, Elisa - Q9Bsi4

Item Number: CM0027415



Introduction

Rabbit polyclonal antibody against human TIN2/TINF2, validated for WB and ELISA. Cross-reacts with mouse. Suitable for shelterin complex and telomere biology research.

[Learn More](#)

Parameter	Value
Product Name	TIN2/TINF2 Rabbit pAb
Remarks/Alias	TIN2; DKCA3; TIN2/TINF2
Species	Human
GeneID (Human)	26277
GeneID	26277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 202-451 of human TIN2/TIN2/TINF2 (NP_001092744.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BSI4
Protein Weight	50kDa
Shipping	Ice bag