

Anti-Tin2/Tinf2 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Bsi4

Item Number: CM0005465



Introduction

Anti-TIN2/TINF2 Rabbit Monoclonal Antibody validated for WB, ELISA; cross-reacts with human, mouse, rat; targets TERF1-interacting nuclear factor 2, a shelterin complex component regulating telomere length and protection.

[Learn More](#)

Parameter	Value
Product Name	TIN2/TINF2 Rabbit mAb
Remarks/Alias	TIN2; DKCA3; TIN2/TINF2
Species	Human
GeneID (Human)	26277
GeneID	26277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-309 of human TIN2/TINF2 (Q9BSI4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BSI4
Protein Weight	50 kDa
Shipping	Ice bag