

Anti-Tiparp Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z3E1

Item Number: CM0016889

Introduction

High-quality rabbit polyclonal antibody targeting human TIPARP (PARP7), validated for WB and ELISA; cross-reacts with mouse. Ideal for studying mono-ADP-ribosylation, AHR signaling, and nuclear protein modifications.

[Learn More](#)

Parameter	Value
Product Name	TIPARP Rabbit pAb
Remarks/Alias	PARP7; ARTD14; pART14; TIPARP
Species	Human
Gene ID (Human)	25976
Gene ID	25976
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 441-657 of human TIPARP (NP_056323.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7Z3E1
Protein Weight	76kDa
Shipping	Ice bag