

Anti-Tir1 Polyclonal Antibody For Wb And Elisa - Q570C0

Item Number: CM0020153



Introduction

Rabbit polyclonal antibody against Arabidopsis TIR1, validated for Western blot and ELISA. Targets the auxin receptor, suitable for plant hormone signaling and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	TIR1 Rabbit pAb
Remarks/Alias	AtTIR1; TRANSPORT INHIBITOR RESPONSE 1; TIR1
Species	Arabidopsis thaliana
GeneID	825473
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 495-594 of arabidopsis thaliana TIR1 (NP_567135.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana
SWISS	Q570C0
Protein Weight	67kDa
Shipping	Ice bag