

Anti-Tirap Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P58753

Item Number: CM0008618



Introduction

High-quality anti-TIRAP rabbit monoclonal antibody for WB, IHC-P, and ELISA. Recognizes human, mouse, rat TIRAP. Explore its role in TLR signaling and innate immune response. Supplied as liquid, ships with ice bag.

[Learn More](#)

Parameter	Value
Product Name	TIRAP Rabbit mAb
Remarks/Alias	Mal, wyatt, BACTS1, MyD88-2, TIRAP
Species	Human
Gene ID (Human)	114609
Gene ID	114609
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TIRAP (P58753).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58753
Protein Weight	24kDa
Shipping	Ice bag