

Anti-Tirap Rabbit Polyclonal Antibody For Ihc-P, Elisa - P58753

Item Number: CM0022833



Introduction

High-quality rabbit polyclonal antibody targeting human TIRAP (Mal), validated for IHC-P and ELISA applications. Ideal for innate immunity, TLR2/4 and RAGE signaling, and inflammatory response studies. SWISS: P58753.

[Learn More](#)

Parameter	Value
Product Name	TIRAP Rabbit pAb
Remarks/Alias	Mal, wyatt, BACTS1, MyD88-2, TIRAP
Species	Human
GeneID (Human)	114609
GeneID	114609
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human TIRAP (NP_001034750.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P58753
Protein Weight	24kDa
Shipping	Ice bag