

Anti-Tirap Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P58753

Item Number: CM0016161



Introduction

Rabbit polyclonal antibody targeting human TIRAP (Mal/MyD88-2), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for innate immunity and TLR/IL-1R signaling pathway studies.

[Learn More](#)

Parameter	Value
Product Name	TIRAP Rabbit pAb
Remarks/Alias	Mal; wyatt; BACTS1; MyD88-2; TIRAP
Species	Human
GeneID (Human)	114609
GeneID	114609
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human TIRAP (NP_001034750.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58753
Protein Weight	24kDa
Shipping	Ice bag