

Anti-Tissue Factor/Coagulation Factor Iii (Cd142) Rabbit Monoclonal Antibody For Wb - P13726

Item Number: CM0006607



Introduction

Rabbit monoclonal antibody against human Tissue Factor (Coagulation factor III, CD142). Validated for WB, IHC-P, IF/ICC, ELISA. Ideal for hemostasis and coagulation research.

[Learn More](#)

Parameter	Value
Product Name	Coagulation factor III/Tissue Factor Rabbit mAb
Remarks/Alias	TF; TFA; CD142; Coagulation factor III/Tissue Factor
Species	Human
Gene ID (Human)	2152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-251 of human Coagulation factor III/Tissue Factor. (NP_001984.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P13726
Protein Weight	33kDa
Shipping	Ice bag