


Anti-Titin (Ttn) Rabbit Polyclonal Antibody For If-P And Elisa - Q8Wz42

Item Number: CM0020511



Introduction

Rabbit polyclonal antibody targeting human Titin (TTN), suitable for IF-P and ELISA, with cross-reactivity to mouse and rat. Useful for studying sarcomere structure, cardiomyopathy, and muscular dystrophy.

[Learn More](#)

Parameter	Value
Product Name	TTN Rabbit pAb
Remarks/Alias	TMD, CMH9, CMD1G, CMPD4, CMYP5, EOMFC, HMERF, MYLK5, SALMY, LGMD2J, LGMDR10, TTN
Species	Human
Gene ID (Human)	7273
Gene ID	7273
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 5347-5604 of human TTN (NP_596870.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WZ42
Protein Weight	3816kDa
Shipping	Ice bag