

Anti-Tkt/Transketolase Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P29401

Item Number: CM0005191



Introduction

Rabbit monoclonal antibody targeting human TKT/Transketolase (P29401), validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying pentose phosphate pathway and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] TKT/Transketolase Rabbit mAb
Remarks/Alias	TK; TKT1; SDDHD; HEL107; HEL-S-48
Species	Human
Gene ID (Human)	7086
Gene ID	7086
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 281-380 of human TKT/Transketolase (NP_001055.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29401
Protein Weight	68kDa
Shipping	Ice bag